

Operating Instructions

LIQUID LIMIT DEVICE

Models: 24-0435 (CL-206)

24-0437 (CL-204A) 24-0441/01 (CL-205-8) 24-0441/02 (CL-205)

ELE International, Inc.

PO Box 608

Loveland, CO 80539

USA

Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Website: www.eleusa.com

Distributor:

ELE International, Ltd.

Chartmoor Road Chartwell Business Park Leighton Buzzard, Bedfordshire

LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249249 E-Mail: ele@eleint.co.uk Website: www.eleint.co.uk

TABLE OF CONTENTS

I. GENERAL INFORMATION	1
II. RELATED USER DOCUMENTATION	1
III. OPERATION	2
IV. SPECIFICATIONS	3
FIGURE 1 MODEL 24-0437 (CL-204)	5
PARTS LIST MODEL 24-0437 (CL-204A)	6
FIGURE 2 MODEL 24-0435 (CL-206)	8
PARTS LIST MODEL 24-0435 (CL-206	9
FIGURE 3 MODEL 24-0441/02 (CL-205)	11
PARTS LIST MODEL 24-0441/02 (CL-205)	12

Rev. 0:10/04

ELE International, Inc. PO Box 608 Loveland, CO 80539 USA Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Distributor: ELE International, Ltd. Chartmoor Road Chartwell Business Park Leighton Buzzard, Bedfordshire LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249200 Fax: + (0) 1525 249249 E-Mail: ele@eleint.co.uk



LIQUID LIMIT DEVICE MODELS 24-0437 (CL-204A), 24-0441/ (CL-205), 24-0435 (CL-206)

I. GENERAL INFORMATION

The Liquid Limit Device, models 24-0437 (CL-204A), 24-0441/ (CL-205) Series, 24-0435 (CL-206), is a device designed to determine the liquid limit of soils. The hard rubber base is formed in a mold to maintain the uniformity of hardness, size and density of all devices. Precision molded cam control parts and nylon bearings provide accurate cup drop and smooth operation. These manufacturing techniques result in a better correlation of test results performed in laboratories around the world.

All models include a removable pin cup assembly for easy cleaning and inspection and ten precision molded **Plastic Grooving Tools** 24-0451 (CL-220). An optional **Grooving Tool** 24-0461 (CL-218A) is available. The 24-0437 (CL-204A) model incorporates a mechanical counter which automatically records the number of cup drops during testing. The 24-0441/ (CL-205) model incorporates a totally enclosed gear motor system which drives the liquid limit device and mechanical counter assembly. The 24-0435 (CL-206) model is a hand operated unit without a mechanical counter.

II. RELATED USER DOCUMENTATION

These operating instructions do not contain the necessary information on the specific test procedures. Please refer to ASTM D-4318 or AASHTO T-89 for additional testing information.

ELE International, Inc. PO Box 608 Loveland, CO 80539 USA Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Website: www.eleusa.com	Distributor:	ELE International, Ltd. Chartmoor Road Chartwell Business Park Leighton Buzzard, Bedfordshire LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249249 E-Mail: ele@eleint.co.uk Website: www.eleint.co.uk
---------------------------------------------------------------------------------------------------------------------------------------------------------	--------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



III. OPERATION

- A. Follow test procedure and operation as given in appropriate documentation.
- B. Position cup to desired distance from base. The hinge position of the cup can be moved forward or backward by loosening the top knurled head screw (Fig. 1, #19; Fig. 2, #13; Fig. 3, #25) and turning the knurled adjusting screw (Fig.1, #1; Fig. 2, #17; Fig. 3, #20). This raises and lowers the cup so its distance from the base is adjustable. After adjustment, secure the desired position by tightening the top knurled head screw.
- C. On models 24-4037 (CL-204A) and 24-0441/ (CL-205) Series, zero the counter by turning the black knob clockwise.
- D. To start the motorized model, 24-0441/ (CL-205) Series, turn the switch to "ON" after connecting the power cord to a power line of appropriate voltage. Model 24-0441/02 (CL-205) requires 110v and Model 24-0411/01 (CL-205-8) requires 220v.
- E. The cup is attached to the sliding holder by a removable pin (Figs. 1 & 2, #3; Fig. 3, #13). Pull the pin out to remove the cup assembly for cleaning and inspection.

ELE International, Inc.
PO Box 608
Loveland, CO 80539
USA
Phone: (970) 663-9780
Fax: (970) 663-9781
E-Mail: soiltest@eleusa.con
Website: www.eleusa.com

Distributor:

ELE International, Ltd.
Chartmoor Road Chartwell Business
Park
Leighton Buzzard, Bedfordshire
LU7 4WG England
Phone: + (0) 1525 249240
Fax: + (0) 1525 249249
E-Mail: ele@eleint.co.uk
Website: www.eleint.co.uk



IV. SPECIFICATIONS

A. MODEL 24-0437 (CL-204A)

Counter Mechanical; records number of cup drops

Cup Spun brass; pin and slot mounting

Base Molded hard rubber

Cam Molded, self-lubricating plastic

Housing Aluminum; holds cup and drop adjustment parts

Grooving Tools ASTM type; plastic; pkg. of 10 (included)

Weight Net 6-1/2 lbs. (2.9 kg)

B. MODEL 24-0441/ (CL-205) SERIES

Gear Motor Attached to metal base; enclosed in housing with on-off

switch

Counter Mechanical; records number of cup drops

Cup Spun brass; pin and slot mounting

Base Molded hard rubber

Cam Molded, self-lubricating plastic

Housing Aluminum; holds cup and drop adjustment parts

Grooving Tools ASTM type; plastic; pkg. of 10 (included)

Distributor: ELE International, Ltd. ELE International, Inc. Chartmoor Road Chartwell Business PO Box 608 Loveland, CO 80539 Leighton Buzzard, Bedfordshire USA LU7 4WG England Phone: (970) 663-9780 Phone: + (0) 1525 249200 Fax: (970) 663-9781 Fax: + (0) 1525 249249 E-Mail: soiltest@eleusa.com E-Mail: ele@eleint.co.uk Website: www.eleusa.com Website: www.eleint.co.uk



Overall Dimensions 8-1/2"w. x 15"d. x 7-1/2"h. (216 x 381 x 191 mm)

Weight Net 14-1/2 lbs. (6.6 kg)

C. MODEL 24-0435 (CL-206)

Cup Spun brass; pin and slot mounting

Base Molded hard rubber

Cam Molded, self-lubricating plastic

Aluminum; holds cup and drop adjustment parts Housing

Hand Crank Aluminum

Grooving Tools ASTM type; plastic; pkg. of 10 (included)

Weight Net 5-1/4 lbs. (3.6 kg)

ELE International, Inc. PO Box 608 Loveland, CO 80539 USA

Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Website: www.eleusa.com

Distributor:

ELE International, Ltd.


Chartmoor Road Chartwell Business

Leighton Buzzard, Bedfordshire LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249249

E-Mail: ele@eleint.co.uk Website: www.eleint.co.uk



FIGURE 1 -- MODEL 24-0437 (CL-204A)



ELE International, Inc. PO Box 608 Loveland, CO 80539 USA Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Website: www.eleusa.com

Distributor:

ELE International, Ltd. Chartmoor Road Chartwell Business Park Leighton Buzzard, Bedfordshire LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249249 E-Mail: ele@eleint.co.uk Website: www.eleint.co.uk



PARTS LIST -- MODEL 24-0437 (CL-204A)

Item No.	Part No.	Description	Req'd.
1	3322-0007	Knurled Adjusting Screw	1
2	3222-0001	Sliding Holder	1
3	9518-0001	4-40 x 3/16 lg. Set Screw	1
4	3346-0001	Crank Wheel Ass'y	1
5	9507-0013	5-40 x 1/4 lg. Soc.Cap Hd. Screw	4
6	1207-0026	Counter	1
7	3346-0003	Hard Rubber Base	1
8	24-0435/10 (CL-227)	Liquid Limit Cup	1
9	3346-0002	Mounting Bracket	1
10	9515-0007	#8 x 1/2 lg. Pan Hd Self-Tap Screw	2
11	9302-0022	Nameplate (Soiltest)	1
12	9401-0019	Rubber Foot	4
13	9302-0035	Data Plate	1
14	9507-0011	1/4-20 x 2 1/4 Flat Head Screw	2

ELE International, Inc. PO Box 608 Loveland, CO 80539 USA Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Website: www.eleusa.com	Distributor:	ELE International, Ltd. Chartmoor Road Chartwell Business Park Leighton Buzzard, Bedfordshire LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249249 E-Mail: ele@eleint.co.uk Website: www.eleint.co.uk
---------------------------------------------------------------------------------------------------------------------------------------------------------	--------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



Item No.	Part No.	Description	Req'd.
15	9302-0028	ELE Logo	1
16	3322-0015	Crankshaft	1
17	9507-0010	5-40 x 1/4 lg. Flat Head Screw	2
18	3322-0003	PIN	1
19	9519-0005	10-32 x 1/2 lg Knurled Screw	1
20	8909-0058	Dowel Pin 3/16 Dia x 1/2 lg.	1
21	3322-0004	Cam Follower	1
22	9517-0028	1/4-20 x 1/4 lg. Set Screw	1
23	3322-0013	Cam	1
24	3322-0014	Housing Ass'y	1
25	24-0451 (CL-220)	Grooving Tool	1 pk.

PO Box 608 Loveland, CO 80539 USA Phone: (970) 663-9780

Phone: (970) 663-9780
Fax: (970) 663-9781
E-Mail: soiltest@eleusa.com
Website: www.eleusa.com


Distributor:

ELE International, Ltd.Chartmoor Road Chartwell Business Park

Park
Leighton Buzzard, Bedfordshire
LU7 4WG England
Phone: + (0) 1525 249200
Fax: + (0) 1525 249249
E-Mail: ele@eleint.co.uk
Website: www.eleint.co.uk



FIGURE 2 -- MODEL 24-0435 (CL-206)



ELE International, Inc. PO Box 608 Loveland, CO 80539 USA Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Website: www.eleusa.com

Distributor:

ELE International, Ltd. Chartmoor Road Chartwell Business Park Leighton Buzzard, Bedfordshire LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249249 E-Mail: ele@eleint.co.uk Website: www.eleint.co.uk



PARTS LIST -- MODEL 24-0435 (CL-206)

Item No.	Part No.	Description	Req'd.
1	3322-0001	Sliding Holder	1
2	3322-0015	Crankshaft	1
3	3322-0003	Pin	1
4	9507-0010	5-40 x 1/4 lg. Flat Head Screw	2
5	3322-0014	Housing Ass'y	1
8	9912-0041	Hard Rubber Base	1
9	9401-0019	Rubber Foot	4
10	3322-0004	Cam Follower	1
11	9517-0025	1/4-20 x 1/4 lg. Set Screw	1
12	24-0435/10	Brass Cup	1
13	9519-0005	10-32 x 1/2 lg Knurled Screw	1
15	3322-0013	Cam	1
16	9507-0011	1/4-20 x 2 1/4 Flat Head Screw	2
17	3322-0007	Knurled Adjusting Screw	1

ELE International, Inc. PO Box 608 Loveland, CO 80539 USA Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Website: www.eleusa.com	Distributor:	ELE International, Ltd. Chartmoor Road Chartwell Business Park Leighton Buzzard, Bedfordshire LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249249 E-Mail: ele@eleint.co.uk Website: www.eleint.co.uk
---------------------------------------------------------------------------------------------------------------------------------------------------------	--------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



Item No.	Part No.	Description	Req'd.
18	9302-0022	Nameplate	1
19	24-0451 (CL-220)	Grooving Tool	1 pk.
20	9302-0035	Data Plate	1

ELE International, Inc.
PO Box 608
Loveland, CO 80539
USA
Phone: (970) 663-9780
Fax: (970) 663-9781
E-Mail: soiltest@eleusa.com
Website: www.eleusa.com Website: www.eleusa.com

Distributor:


ELE International, Ltd.


Chartmoor Road Chartwell Business Park

Park
Leighton Buzzard, Bedfordshire
LU7 4WG England
Phone: + (0) 1525 249200
Fax: + (0) 1525 249249
E-Mail: ele@eleint.co.uk
Website: www.eleint.co.uk



FIGURE 3 -- MODEL 24-0441/02 (CL-205)





ELE International, Inc. PO Box 608 Loveland, CO 80539 USA Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Website: www.eleusa.com Pistributor: ELE International, Ltd. Chartmoor Road Chartwell Business Park Leighton Buzzard, Bedfordshire LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249249 E-Mail: ele@eleint.co.uk Website: www.eleint.co.uk



PARTS LIST -- MODEL 24-0411/02 (CL-205)

Item No.	Part No.	Description	Req'd.
1	3302-0006	Shaft	1
2	1207-0014	Counter	1
3	3302-0007	Cover	1
4	3302-0010	Mounting Base	1
5	8805-0003	Motor	1
6	3302-0018	Sprocket, 10 Teeth 5/16 Bore	1
7	3302-0011	Motor Mounting Bracket	1
8	9427-0008	Chain 1/4 Pitch, 33 Link	1
9	9427-0012	Chain Master Link, 1/4 Pitch	1
10	3346-0002	Counter Mounting Bracket	1
11	3302-0017	Sprocket 12 Teeth 1/4 Bore 24-0441/01 (CL-205-8)	1
12	3322-0004	Cam Follower	1
13	3322-0003	Pin	1
14	3322-0001	Sliding Holder	1
15	0907-0027	Toggle Switch	1

ELE International, Inc. PO Box 608 Loveland, CO 80539 USA Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Website: www.eleusa.com	ELE International, Ltd. Chartmoor Road Chartwell Business Park Leighton Buzzard, Bedfordshire LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249249 E-Mail: ele@eleint.co.uk Website: www.eleint.co.uk
---------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



16	9403-0003	Line Cord Grommet	1
17	3302-0022	Housing Ass'y.	1
18	24-0435/10 (CL-206)	Brass Cup	1
19	3302-0012	Hard Rubber Base	1
20	3322-0007	Knurled Adjusting Screw	1
21	9517-0028	1/4-20 x 1/4 lg. Set Screw	1
22	3322-0013	Cam	1
23	9507-0010	5-40 x 1/4 lg. Flat Head Screw	2
24	9507-0011	1/4-20 x 2 1/4 Flat Head Screw	2
25	9519-0005	10-32 x 1/2 lg Knurled Screw	1
27	9507-0013	5-40 x 1/4 lg. Soc. Cap Hd. Screw	4
28	9302-0022	Nameplate	1
29	3302-0015	Sprocket 1/4 Bore 10 Teeth	1
30	9401-0019	Rubber Feet	4
31	9507-0012	1/4-20 x 1 lg. Flat Hd. Mach. Screw	2
32	8903-0003	Counter Coupling	1

ELE International, Inc. PO Box 608 Loveland, CO 80539 USA Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Website: www.eleusa.com	Distributor:	ELE International, Ltd. Chartmoor Road Chartwell Business Park Leighton Buzzard, Bedfordshire LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249249 E-Mail: ele@eleint.co.uk Website: www.eleint.co.uk
---------------------------------------------------------------------------------------------------------------------------------------------------------	--------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



Item No.	Part No.	Description	Req'd.
33	3302-0013	Shaft Spacer	2
35	9515-0007	#8-1/2 lg. Pan Hd. Type A Philips, Self Tapping Screw	2
36	8608-0015 8608-0016	3-Conductor Line Cord w/Plug 24-0441/02(CL-205) 24-0441/01 (CL-205-8)	1
37	2205-0003	Wire Nut	3
38	9506-0024	#10-32 x 3/4 lg. Flt. Hd. Mach. Screw	2
39	9505-0002	#8-32 x 1/2 lg. Soc. Hd. Cap Screw	2
40	0907-0001	On-Off Switch Plate	1
41	9411-0002	Cable Clamp	1
42	9511-0012	#8-32 x 1/2 lg. Soc. Hd. Cap Screw	2
43	9701-0018	#8 Brass Washer	2
44	9702-0003	#8 Lock Washer	2
45	1901-0064	"Push On" Terminal	4
46	8601-1304	Stranded Wire #20 x 2 lg. Blue	1
47	8601-1004	Stranded Wire #20 x 6 lg. Black	1
48	8702-0007	#18 Plastic/Glass Tubing 4" lg.	1

ELE International, Inc. PO Box 608 Loveland, CO 80539 USA Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Website: www.eleusa.com	Distributor:	ELE International, Ltd. Chartmoor Road Chartwell Business Park Leighton Buzzard, Bedfordshire LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249249 E-Mail: ele@eleint.co.uk Website: www.eleint.co.uk
---------------------------------------------------------------------------------------------------------------------------------------------------------	--------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



Item No.	Part No.	Description	Req'd.
49	1508-0005	Transformer 24-0441/01 (CL-205-8)	1
50	8601-0029	Wire #18 Black 1' 24-0441/01 (CL-205-8)	1
51	24-0451 (CL-220)	Grooving Tool	1 pk.
52	9518-0023	#10-32 x 3/16 Soc Set Screw	2
53	1901-0034	Terminal #6	1
54	9302-0035	Data Plate	1

ELE International, Inc.
PO Box 608
Loveland, CO 80539
USA
Phone: (070) 663 0790

Phone: (970) 663-9780 Fax: (970) 663-9781 E-Mail: soiltest@eleusa.com Website: www.eleusa.com Distributor:

ELE International, Ltd.
Chartmoor Road Chartwell Business

Park Leighton Buzzard, Bedfordshire LU7 4WG England Phone: + (0) 1525 249200 Fax: + (0) 1525 249249

Fax: + (0) 1525 249249 E-Mail: ele@eleint.co.uk Website: www.eleint.co.uk